

**Tri-Power Enzymatic Cleaner**  
**Instructions for Use**  
**ASP100**

Tri-Power Enzymatic Cleaner is a true all-in-one product which has been approved for use in every step of the decontamination process. This includes the sink (hand-washing), ultrasonic, and for all cycles in the washers/sterilizers. The only additional product that may be needed, depending on hospital protocol, is an instrument lubricant at the end of the decontamination process. Tri-Power is completely free-rinsing and no residue will adhere to the instruments after use. It is very important that Tri-Power is not mixed with other chemicals/products. Tri-Power Enzymatic Cleaner's manufacturing process is ISO9001:2000 certified and is it registered with the FDA as a Class 1 medical device.



**Instructions for Use**

**Manual-wash**

Soak time: 3 to 5 minutes.

Dilution rate: 1/8 to 1oz of solution per gallon of water.

Rinse thoroughly with water.

May be used in all water temperatures but is recommended to be used in as warm water as tolerable by the person hand-washing the instruments. (There is no need to monitor water temperature)

**Ultrasonic Cleaner**

Dilution rate: 1/8 to 1oz of solution per gallon of water.

Rinse thoroughly with water.

May be used in all water temperatures. Optimal temperature is in the range of 70-140 degrees. Change solution when visibly soiled.

**Washer/Disinfector – Standard Concentration**

Dilution rate: 1/8 to 1oz of solution per gallon of water.

May be used in all water temperatures. The optimal temperature is 140 degrees.

Cold water should not be used during the wash cycle.

**Cautions / Warnings / Stability**

Keep out of reach of children. Avoid prolonged contact with skin. Irritant to eyes and skin. Wear gloves when handling. In case of eye contact, immediately flush with water and consult with a physician, if necessary. Do not ingest. Discard used solution. Do not reuse. Tri-Power is not a sterilizer or disinfectant. Instruments must be followed up with standard sterilization procedures. Refer to SDS for safety recommendations. Store in a cool dry place. **18 month Product Stability at or below room temperature and uncontaminated.**