

Tri-Power Enzymatic Cart Wash

Tri-Power Enzymatic Cart Wash is completely free-rinsing and no residue will adhere after use. It is non-corrosive, suitable for use on all materials, eliminates rusting, scaling, spotting and keeps carts and trays looking new. It is very important that Tri-Power is not mixed with other chemicals/products. Tri-Power Enzymatic Cleaner's manufacturing process is ISO9001:2000 certified and is it registered with the FDA as a Class 1 medical device.

Instructions for Use

Cart Washer:

Follow manufacturer's recommended instructions

Dilution rate: 1/2 to 1 oz of solution per gallon of water.

Can be manually or automatically dispensed into cart washer