

## **Tri-Power Rinse Aid**

Tri-Power Rinse Aid is a concentrated drying aid for use in the final rinse of cart washers which speeds the drying process. Tri-Power Rinse Aid promotes the free rinsing of water to resist spotting on surfaces which may cause metal surface corrosion. This product is designed for use on glassware, surgical instruments, rubber, metal and plastic items in automated washers only. It is non-corrosive, suitable for use on all materials, eliminates rusting, scaling, spotting and keeps carts and trays looking new. The Tri-Power Rinse Aid manufacturing process is ISO9001:2000 certified and is it registered with the FDA as a Class 1 medical device.

### Instructions for Use

#### **Cart Washer:**

Follow manufacturer's recommended instructions

Dilution rate: 1 to 3oz of solution per gallon of water during rinse cycle.

Can be manually or automatically dispensed into cart washer

